



Improve your quality of life after medially stabilized knee arthroplasty





IMPROVE YOUR QUALI



TY OF LIFE

	Introduction	5
1	The knee and the arthritis Discover how your knee works and how arthritis may affect you	6
2	Getting ready for your operation Things to do to prepare for your operation	8
3	Total knee replacement Learn about an artificial knee and how it may benefit you	10
4	GMK Sphere: stability for life Discover why GMK Sphere total knee replacement is different	12
5	While you are in the hospital Your stay in hospital until you are discharged	14
6	Taking care of your new kneeenjoy your new knee!	16

This brochure has been produced to help you feel safe and confident about your operation. Questions you may have regarding the surgery and post-operative recovery, are addressed.



IMPROVE YOUR QUALI



TY OF LIFE

INTRODUCTION

The knee is the largest and most complex joint of our body. It has a very difficult job; carrying the weight of our body with every step we take! Therefore, it is not surprising that knees are the joints in our body that are most vulnerable to injuries or developing degenerative joint diseases, such as arthritis.

One of the consequences of any joint disease is pain.

Statistics show that roughly one third of the American population over the age of 45 suffers from knee pain. Knee pain may limit your daily activities, affect your fitness level, emotional health, and your **general well-being!**

Can you get rid of the pain?

There are a variety of surgical and non surgical solutions to treat your disease. Your doctor will prescribe the most suitable treatment, according to your age, activity level and expectations.

Knee pain and stiffness caused by advanced arthritis is severely limiting and your physician may suggest you undergo a total knee replacement.




1 - THE KNEE AND TH





E ARTHRITIS

ANATOMY OF THE KNEE

The knee joint is comprised of three bones: the thigh bone (femur), the shin bone (tibia) and the kneecap (patella). The leg movement is driven by the thigh muscles, the biggest one being the quadriceps, located in the front of the thigh. The thigh and shin bones are connected by ligaments, which give stability to the knee joint. The surface of the kneecap, thigh bone and shin bone, where the bones come in contact, is coated with a smooth tissue called cartilage. The cartilage, together with a substance called synovial fluid, prevents the bones from rubbing against each other and causing damage.









ARTHRITIS OF THE KNEE

With arthritis, the cartilage deteriorates and the bones start rubbing directly against each other. The result is joint **pain**, which can become worse over time, and limits motion.

Knee replacement is a common treatment for severe arthritis.

Successful knee replacement surgery can result in dramatic pain relief and improvement in knee joint function.







2 - GETTING READY FO



MEDICAL EXAM

Before undergoing your total knee replacement, your doctor will prescribe a complete physical examination to assess your condition and to evaluate factors that could potentially interfere with your surgery.

r your operation

TESTS

Your doctor may ask you for blood and urine tests and possibly a cardiogram, prior to your surgery.

CHECK YOURSELF

Check your skin

If your knee and leg have any skin infections or irritation, contact your orthopaedic surgeon prior to surgery: he or she will tell you how best to prepare your skin for surgery.

2 Check your teeth

The incidence of infection after knee replacement is very low, but even an infection far from the knee joint can spread through the bloodstream to the new knee. Prevention is the best way to avoid early problems. Therefore, you should contact your dentist to have your teeth checked before your Total Knee Arthroplasty surgery.

MEDICATIONS

Prior to surgery, provide your surgeon with a complete list of the medications you are taking including doses and times. He or she will inform you if you need to stop or change any medication.

SPECIAL EQUIPMENT

After the surgery, special equipment, such as support stockings and crutches, may be needed. Your doctor and therapist can assist you with obtaining these aids.



3-TOTAL KNEE REPLA

WHAT IS TOTAL KNEE REPLACEMENT?

Total knee replacement surgery aims to substitute the bone and cartilage of the joint damaged by arthritis with metallic and plastic implants.

The surfaces of the thigh and shin bones are replaced with high-resistant metallic components, called the **femoral component** and **tibial baseplate**.

Between the femoral component and the tibial baseplate, a **plastic tibial insert** is implanted. It replaces the cartilage function allowing the thigh and shin bone to slide on each other. All materials used in a total knee replacement are highly biocompatible.







PLASTIC TIBIAL INSERT

Polyethylene



TIBIAL BASEPLATE

Metal



CEMENT

WHY TOTAL KNEE REPLACEMENT?

With almost 50 years of history, total knee replacement surgery is a very common and safe procedure for the treatment of severe arthritis. Approximately 1,000,000 knee replacements are performed annually worldwide. The main benefits of a successful total knee replacement are:

Reduction of knee pain
Pain may be rapidly and

dramatically reduced, potentially eliminated!

Recovery of mobility
You may dramatically improve the mobility of your knee.

3 Improvement in quality of life

Your everyday activities may no longer be limited by pain and reduced mobility!






4-GMK SPHERE: STAB

GMK SPHERE: THE MEDIALLY STABILIZED KNEE

A normal knee is a complex structure: when it bends, the lateral (outer) side rolls back while the medial (inner) side remains stable.

The GMK Sphere knee is designed to provide a more natural motion that replicates the movement of the healthy knee: it allows **stability** in the **medial side** through a ball-in-socket mechanism and accommodates the **natural movement** of the **lateral side** without imposing any constrained motion.



STABILITY IS IMPORTANT FOR PATIENTS


The GMK Sphere medial ball-in-socket mechanism can replicate the function and the stability provided by the natural stabilizers of the knee, most of which are removed during total knee replacement.

Stability is important for common activities such as ascending and descending stairs, standing up from a chair, getting into a car as well as more demanding activities like shopping and gardening.

Published studies* actually confirm that **patients prefer knees with a medial ball-in-socket** compared to other conventional knee designs, as they feel more natural, more stable and stronger during daily activities.

^{*} Pritchett JW, "Patients Prefer A Bicruciate-Retaining or the Medial Pivot Total Knee Prosthesis", J. Arth 2011.

LITY FOR LIFE





5 - WHILE YOU ARE IN



TAKE WITH YOU

THE HOSPITAL

THE SURGICAL PROCEDURE

When the surgical team is ready, you'll be taken to the operating room. The surgery may be done through a mid-vastus (through the front) or a vastus (off to the side) approach to the knee joint. The damaged bone is removed from the articulating surfaces of both the femur (thigh bone) and tibia (shin bone) so that the new metal prosthesis can be fitted. After the new joint is in place, skin incisions are closed with staples, stitches, or surgical glue.

THE DAY OF YOUR OPERATION

The surgical procedure usually takes 1 to 2 hours.

Surgery will begin with pre-surgical preparation and followed by monitoring in the recovery room. The time away from your room will be longer than the operation due to time needed for your preparation for surgery, administration of anaesthesia and monitoring as you recover from the anaesthesia. Careful measures will be taken to relieve pain that may occur after the surgery. Regular checks will be made by the nurses.

AFTER THE OPERATION

Doctors, nurses and therapists will take care of your recovery by defining the most suitable rehabilitation program for you and accompanying you through the gradual recovery process.

Rehabilitation may begin the day of the operation, subject to your doctor's approval. However, you will progress to weight bearing activities as tolerated and may discontinue assistive devices as your comfort level improves, always following your doctor's prescriptions.

- A complete list of your routine medications including doses, times and your x-rays
- 2 All papers for hospital admission including a copy of insurance cards
- 3 Any equipment (crutches, stockings, etc.) ordered for you by your doctor



6-TAKING CARE OF Y

LONG TERM CARE OF YOUR NEW KNEE

Follow your orthopedic surgeon's instructions carefully to minimize any potential complications which may affect your recovery. These complications are quite infrequent and some simple rules can dramatically reduce their likelihood.

DON'T FORGET

- Lead a healthy and active life.
- 2 In case of fever, throat inflammation, pulmonary inflammation or similar, tell your physician that you have a knee implant.
- 3 Undergo regular general check-ups.

OUR NEW KNEE

If you have any concerns about your new knee, don't hesitate to contact your physician and, finally, ...

...enjoy your new knee!






IMPROVE YOUR QUALI

My experience:		



TY OF LIFE



For further information visit the website:

medacta.com

