


SecureFix™

PEEK ALL-INSIDE MENISCAL REPAIR

SecureFix™ XL

PEEK ALL-INSIDE MENISCAL REPAIR



UHMWPE SUTURE PEEK BUTTON MENISCUS REPAIR IMPLANT

**DESIGNED FOR ACCURACY, RELIABILITY
AND EASE OF USE IN SURGICAL PROCEDURES**

ERGONOMIC DESIGN


enables one-handed use
with a 3-step deployment process

AVAILABLE IN MULTIPLE CONFIGURATIONS

Pre-bent Up, Down, Left, Right or Flexible
(up to 45 degrees) needle*

*Available also in XL needle configuration


UHMWPE SUTURE PEEK BUTTON MENISCUS REPAIR IMPLANT



3-STEP DEPLOYMENT PROCESS



PIERCE
the needle



TURN
the safety knob



FIRE
the implant

ORDERING INFORMATION

Order Number	Description	Order Number	Description
05.06.0010	SecureFix Meniscus Repair Implant	05.06.0020	SecureFix XL Meniscus Repair Implant
05.06.0011	SecureFix Meniscus Repair Implant Left	05.06.0021	SecureFix XL Meniscus Repair Implant Left
05.06.0012	SecureFix Meniscus Repair Implant Right	05.06.0022	SecureFix XL Meniscus Repair Implant Right
05.06.0013	SecureFix Meniscus Repair Implant Down	05.06.0023	SecureFix XL Meniscus Repair Implant Down
05.06.0014	SecureFix Meniscus Repair Implant Flexi	05.06.0024	SecureFix XL Meniscus Repair Implant Flexi

ASSOCIATED INSTRUMENTS

Order Number	Description	Disposable/Reusable
05.06.10.0110	SecureFix - Meniscal open cannula(sterile)	Disposable
05.06.10.0111	SecureFix -Meniscal knot pusher cutter(sterile)	Disposable
05.06.10.0100	SecureFix - Meniscal open cannula	Reusable
05.06.10.0101	SecureFix - Meniscal knot pusher cutter	Reusable
05.06.10.0102	SecureFix - Meniscus Diamond Rasp	Reusable
05.06.10.0103	SecureFix - Meniscus Depth Gauge	Reusable

