


MedactaLock PEEK RC

J-LOCK PEEK SUTURE ANCHOR



KNOTLESS ANCHOR FOR ROTATOR CUFF REPAIR

NO DRILLING REQUIRED


Pilot hole is created with a punch
Anchor is malleted in afterwards

TECHNICAL DATA

IMPLANT LENGTH: 15mm
Available in 3 sizes:
4.5mm/ 5.5mm/ 6.5mm

**SUITABLE FOR LATERAL ROW
ROTATOR CUFF REPAIR**

**SUTURE EYELET
INTEGRATED IN ANCHOR**



ORDERING INFORMATION

MECTALOCK PEEK RC			
Order Number	Description	Order Number	Description
05.10.004	MectaLock PEEK – Ø4.5 - Short driver	05.10.006	MectaLock PEEK – Ø6.5 - Short driver
05.10.005	MectaLock PEEK – Ø5.5 - Short driver		

DEDICATED INSTRUMENTS			
Order Number	Description	Order Number	Description
05.10.10.0021	Punch – Ø4.5	05.10.10.0023	Punch – Ø6.5
05.10.10.0022	Punch – Ø5.5		

RELATED PRODUCTS

Order Number	Description	Order Number	Description
05.15.10.0015	FastShuttle - Suture Passer & Retriever	05.14.10.0001	Shoulder Cannula System Ø8.5x70mm
05.15.10.0016	FastShuttle - Suture Passer	05.14.10.0002	Shoulder Cannula System Ø8.5x90mm
05.15.10.0017	FastShuttle - Disposable Needle - STERILE (5ea)	05.14.10.0035	Shoulder Cannula System Ø8.5x70mm for Ø4mm Switching Stick
05.10S.001	Shoulder General Set	05.14.10.0036	Shoulder Cannula System Ø8.5x90mm for Ø4mm Switching Stick

